

INDEX

INDEX	1
SYMBOLS AND ABBREVIATIONS	5
INTRODUCTION	7
pH homeostasis	7
pH regulation mechanism	8
<u>Intracellular buffering</u>	<u>8</u>
<u>Biochemical pH-stat</u>	<u>9</u>
<u>Biophysical pH-stat</u>	<u>10</u>
<u>H⁺-ATPases</u>	<u>10</u>
<u>ATPases and secondary active transport</u>	<u>12</u>
Regulation of pH homeostasis in <i>Saccharomyces cerevisiae</i>	13
Weak acids and pH homeostasis	13
Mechanisms of weak acid toxicity in microorganisms	14
Acid stress and pH signaling in yeast	16
The soil acidity	17
pH signaling in <i>Arabidopsis thaliana</i>	17
<u>Plant cell growth</u>	<u>17</u>
K⁺ transport in <i>Arabidopsis thaliana</i>	18
K⁺ channels	19
<u>Shaker channels</u>	<u>19</u>
<u>Inwardly Rectifying Conductance (IRC)</u>	<u>21</u>
<u>Weakly Inwardly Rectifying Conductance (WIRC)</u>	<u>22</u>
<u>Outwardly Rectifying Conductance (ORC)</u>	<u>22</u>
<u>KCO-2P channels (K⁺ Channel Ca²⁺ Outward rectifying 2 Pore)</u>	<u>23</u>
<u>Kir-like channels (K⁺ inward rectifying like)</u>	<u>23</u>
<u>Non-voltage-gated channels</u>	<u>23</u>
<u>GLR (Glutamate Receptors)</u>	<u>23</u>
<u>CNGC Channels (Cyclic Nucleotide-Gated Channel)</u>	<u>23</u>
K⁺ transporters	24
<u>HAK (High Affinity K⁺ Uptake)/ KUP (K⁺ Uptake Permease)</u>	<u>24</u>
<u>Antiporter family</u>	<u>25</u>
<u>Family of high affinity K⁺ transporters</u>	<u>26</u>

H⁺-ATPase and K⁺ transport can play more roles	26
<u>Mislocalization of ABCB1, ABCB4, ABCB19</u>	<u>28</u>
Auxins	29
<u>IAA synthesis</u>	<u>29</u>
<u>Long distance IAA movements</u>	<u>29</u>
<u>Cell to cell movement</u>	<u>31</u>
<u>IAA genetic response</u>	<u>32</u>
Different auxins, differen roles	33
Immunophilins	35
FKBPs in yeast and <i>Arabidopsis thaliana</i>	36
<u>Enzymatic activity and its substrates</u>	<u>37</u>
ROF1 (FKBP62) and ROF2 (FKBP65)	41
Heat Shock Proteins and co-chaperons	42
<i>Arabidopsis thaliana</i>, <i>Saccharomyces cerevisiae</i>: two model systems	44
OBJECTIVES	47
MATERIAL AND METHODS	49
Plasmid maps	49
<i>Arabidopsis</i> solutions	50
<u>Murashige-Skooge medium (MS)</u>	<u>50</u>
<u>K⁺-free medium</u>	<u>50</u>
<i>Arabidopsis</i> methods	51
<u>Plant material, growth and transformations</u>	<u>51</u>
<u>Selection of At5g48570 and At3g25230 knock-out <i>Arabidopsis</i> lines</u>	<u>51</u>
<u>PCR reaction (genotyping)</u>	<u>51</u>
<u>Over-expression of ROF2 in transgenic plants</u>	<u>52</u>
<u>K⁺ Determination</u>	<u>52</u>
<u>Measurement of Rb⁺ uptake</u>	<u>53</u>
<u>Measurement of PM H⁺ ATPase activity in <i>Arabidopsis</i> roots</u>	<u>53</u>
<u>Root acidification assay</u>	<u>53</u>
<u>Electrophysiology of <i>Arabidopsis</i> root epidermal cells</u>	<u>54</u>
<u>Fluorescence microscopy of ROF2 GFP fusion protein</u>	<u>55</u>
<u>Hormonal determination</u>	<u>55</u>
RESULTS	57

Initial data	57
<i>in silico</i> analysis	57
Plant immunophilin mutant	61
<u>Loss-of-function</u>	61
<u>Gain-of-function</u>	62
Positive control in germination assay	63
Germination in weak acid stress conditions	65
Growth in presence of acetic acid in liquid medium	67
ATPase activity	68
Tolerance with toxic cations	69
Germination in low K⁺ medium	72
K⁺ uptake	73
H⁺ extrusion	74
Electrophysiological measurements	75
Sub-cellular localization during intracellular acid stress	77
Flowering delay	78
Apical dominance and seeds	79
Root branching and hormonal level	80
Cold sensitivity	82
Post-transcriptional regulation	83
DISCUSSION	87
CONCLUSIONS	93
BIBLIOGRAPHY	95
APPENDIX	107
Rutinary protocols	107
<u><i>E. coli</i> solutions</u>	107
<u>Antibiotics</u>	107
<u>LB (Luria-Bertani) medium</u>	107
<u>X-gal: 5-Bromo-4-chloro-3-indolyl-β-D-galactoside</u>	107
<u>IPTG: Isopropylthio-β-D-galactoside</u>	107
<u><i>E. coli</i> methods</u>	107

<u>Plasmid purification (miniprep)</u>	<u>107</u>
<u>Second DNA purification (electrophoration grade)</u>	<u>108</u>
<u>Competent cells (Heat-Shock)</u>	<u>108</u>
<u>Transformation (Heat-Shock)</u>	<u>108</u>
<u>Competent cells (Electrophoration)</u>	<u>109</u>
<u>Transformation (Electrophoration)</u>	<u>109</u>
<u><i>Agrobacterium tumefaciens</i> methods</u>	<u>109</u>
<u>Competent cells (Electrophoration)</u>	<u>109</u>
<u>Transformation (Electrophoration)</u>	<u>110</u>
<u><i>Arabidopsis</i> protocols</u>	<u>110</u>
<u>Seed sterilization</u>	<u>110</u>
<u>Plant RNA extraction</u>	<u>110</u>
<u><i>Arabidopsis</i> transformation by floral dip</u>	<u>110</u>
<u>Genomic DNA Extraction with CTAB (Soltis' Lab Protocol)</u>	<u>111</u>
<u><i>in vivo</i> ATPase activity assay</u>	<u>112</u>
<u>RNA and DNA Methods</u>	<u>113</u>
<u>DNA Electrophoresis</u>	<u>113</u>
<u>RNA Electrophoresis</u>	<u>113</u>
<u>Denaturant gel</u>	<u>113</u>
<u>No denaturant gel</u>	<u>114</u>
<u>Transferring RNA to membrane</u>	<u>114</u>
<u>Radiolabelling probe</u>	<u>114</u>
<u>Transferring northern blot pre-hybridization and hybridation</u>	<u>115</u>